

## Powder sampler Mini



The Powder Sampler Mini is a stainless steel sampler for small quantities.

When the button is pressed, a sampling rod enters the sample material and retracts when released. In this way, a fixed amount of sample material is enclosed in the sample notch of the rod and removed.

The sampler precisely collects 0.01 ml, 0.05 ml and 0.1 ml sample volumes, depending on the version.

With an insertion depth of up to 300 mm, multiple target point samples can be collected at different depths.

The small-volume powder sampler is made of high-quality stainless steel 1.4404/1.4571 (316L/316Ti).

- ▶ Sampling of small quantities (0.01 - 0.1 ml)
- ▶ Simple and fast application with the push of a button

Powder sampler Mini					
Chamber volume ml	Length mm	Dia. mm	Item no.		
0,01	300	6	5346-3010	N4	
0,05	300	6	5346-3050	N4	
0,1	300	6	5346-3100	N4	

## PharmaSpoon | for ultrapure sampling



PharmaSpoon for taking samples of powdery substances, intermediate products and material in the pharmaceutical industry. Special materials with excellent surface qualities must be used for this to prevent cross-contamination. Our PharmaSpoon fulfils these requirements. Made of stainless steel AISI 316L (1.4404), surface highly polished by hand (not electro-polished!). The excellent surface properties of the PharmaSpoon provide significant advantages.

- ▶ Possible residues are even visible in the trace sector ( $\mu\text{g}$  sector)
- ▶ The surface is easy to clean
- ▶ The PharmaSpoon is also suitable for sterile sampling

The PharmaSpoon was developed by an internationally known pharmaceutical company in Basel, Switzerland, and its usefulness has been confirmed in actual practice.

PharmaSpoon					
Total length mm	Item no.				
650	5324-0001				N2